leonardo

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***System created
to control
movement
and organize
therapeutic exercises***

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**System created
to control
movement and
organize therapeutic
exercises**

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The term motor learning defines the ability to learn or modify motor skills in a relatively permanent way, by means of a complex perceptual - motor - cognitive process that leads to the individual's increased performance as a result of their practice and interaction with the environment (Magill, 1980; Woollacott & Shumway-Cook, 1994; Schmidt & Lee, 2005). Rehabilitation is actually a true motor re-learning process (Perfetti, 1979; Carr & Shepherd, 1987) aimed at acquiring new skills and recovering motor capabilities that were lost by the individual suffering from any conditions. An essential condition for learning is the existence of feedback to return intrinsic or extrinsic information on the outcome and result from one's motor response (Adams, 1987; Winstein, 1991). This proves to be true both upon developing a movement pattern and when matching the motion produced with its correct reference. The performance we require from the patient must be gradual and suited to their capabilities, so as to generate continual improvement and progressive adjustment in the physician's requests. Depending on the source originating the sensory feedback, we can talk of intrinsic feedback and augmented feedback or extrinsic feedback (Winstein, 1987). Extrinsic feedback is given by information coming from an external source to the individual who is performing the movement and that information, in turn, increases and amplifies the intrinsic feedback. This type of information is always linked to how movement is executed or to particular aspects of it that the individual cannot or is unable to be aware of. It can be provided through various types of channels having different sources and produced in different variations. It is now widely recognized that such information has a key role in providing acknowledgment as to whether the movement is correct and reinforcing learning as a result (Schmidt & Lee, 2005; Wulf, 2007; Roelfsema et al., 2010). It is in this context that I have developed a series of tools over the years that target the patient's cognitive processes, particularly perception, attention and motor movement control (Crippa, 2009), both with regard to central or peripheral damages. It is important to guide the patient towards an advanced re-learning of gestures: it is not about a purposeless and rough movement, but about searching for recruitment quality while controlling compensation and elementary schemes.

The following exercises certainly do not purport to be a path to be followed to set up a treatment program, but they rather represent a number of stimuli that, before implementation, require personal education and a careful evaluation of the patient in order to define a correct rehabilitation path.

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LEONARDO MODULAR BOARD

This aid can be used in a rehabilitation program to gradually restore proprioception through various exercises.

The board can be attached to hinges, half-spheres or spherical fittings through a series of holes and magnetized pins, which makes it possible to use it for a variety of exercises.

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LEONARDO POSITIONAL FEEDBACK

These positional feedback systems are especially useful as assessment and rehabilitation tools.

They allow to monitor and amplify movements (of the head, trunk, pelvis and knee) in order to provide a quick assessment of asymmetries, compensations and changes in patient performance. Once the problem has been identified, it is possible to target the action on the chosen body segment, guiding the patient towards the acquisition of new strategies.

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LEONARDO COMPLETE

This kit includes the same accessories found in the Leonardo Modular Board and Positional Feedback systems, plus another series of valuable aids for the rehabilitation process in orthopedic and neurological pathologies that help design specific exercise program for the patient's pathology.

The kit has been designed to control movement and plan therapy exercises, and provides an immediate answer to the various therapeutic needs of the daily gym exercise.

Thanks to its modular structure, it allows to create numerous therapeutic aids, significantly reducing the quantity of tools required. This ensures an opportunity for professional intervention even at the patient's home.

The special connecting systems supplied make it possible to easily set up the various therapeutic aids.

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The physiotherapist asks the patient to identify or assume different positions by adducting/abducting his/her upper limb, which is placed on the board. Thanks to bearing spheres, the task is significantly facilitated, allowing even the smallest recruitment to be indicated.

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Initially, the primary task the patient is asked to perform is to let him/herself be guided by the physiotherapist while identifying the different positions and with his/her hand placed on the board. In this manner, the upper limb is used as an extension of the trunk. Afterwards, the patient will actively assume different positions. The inclinometer placed on the iliac crest controls pelvis rotation. In this manner, it is possible to create tasks to release the girdle joints.

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The patient, who is in a supine decubitus position, places the upper portion of his/her trunk and his/her pelvis on the 2 oscillating boards. While keeping the proximal board stable, the patient is asked to reach the different levels with the board under his/her pelvis. The different levels are graphically indicated on a graduated scale. Failure for the trunk to fragment activates acoustic FB.

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The patient is asked for the refined recruitment of his/her knee extensors to identify the weight located at the extremity opposite the oscillating board.

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The forearm is placed in an intermediate position. The patient is directed to follow the arch of the circumference that refers to the thumb contrarotation.

2 snap-on pins allow fastening the aid to the board and selecting among the 2 trajectories.

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By drawing the 2 boards closer and inserting the panel pins in the appropriate holes, the patient is asked to identify the positions of the forearm in a prone-supination position.

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The patient's hand is placed on the oscillating board. The different holes both on the 2 sides and at the base of the surface allow personalizing the position of the rotating axis. The pointer on the board allows separating flex/extension movements by the wrist and identifying them with reference to height.

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
The patient is asked to identify the different shims located under the proximal edge of the board. The pyrographed number on each cylinder indicates height.

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Exercise similar to number 7. The simple pointer has been replaced by a laser pointer. The latter, placed on the oscillating board, allows magnifying wrist flexion/extension movements.

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By learning to control finger distal irradiation, the patient increases pressure on his/her palm to identify the resistance by the number of springs located under the oscillating board.

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The patient is asked to use his/her upper limb to perform the requested trajectory by moving the board on the work surface while maintaining or achieving the requested wrist position. The extent of recruitment depends on the use of springs.

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Fingers' weighing request while the patient's hand is placed on the small oscillating board. The aid's holes allow identifying the appropriate height and lever arm. The physiotherapist inserts and removes washers, which the patient is asked to quantify.

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In this exercise, with the insertion of the thin indicator, the main task consists in achieving/identifying the flexion-extension finger positions. The ability to modify the weight by adding washers ensures exercise variability.

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The patient's hand is placed on the board with a central pin and using the pointer, follows the trajectory determined by the physiotherapist. A metal contact connected both to the pointer and the flexible wire closes the circuit while detecting any errors or inaccuracies in performing the task.

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Abduction of the upper limb with the hand resting on the board positioned on the magnetized spheres. The thoracic indicator connected to the FB system, with its optimal calibration range, signals the failed dissociation between the trunk and arm

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Abduction of the upper limb with the elbow extended on the horizontal plane performing kinaesthetic control tasks with limb weighing request. The thoracic indicator connected to the FB system shows failed dissociation between the trunk and the arm.

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This aid consists of a board equipped with bearing spheres, on which a second board is positioned, whose holes are facing upward. A hemisphere is inserted in one of the holes. The hand of the upper, plegic limb is placed on the aid. The physiotherapist moves the latter onto the work surface. The patient must reach the hemisphere using his/her healthy hand. This work allows gathering and processing information related to the metric space and somesthetic relations between limbs

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The level connected to the oscillating board controls the horizontal position of the board by continuously making tonic adjustments and performing a balancing action between the two limbs

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
The aid consists of a board equipped with bearing spheres, over which a second board is positioned with its transversal axis connected to the inclinometer. The patient is asked to keep the horizontal position of the oscillating board while extending the elbows and moving the trunk forward.

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The laser applied to the forearm magnifies elbow angle excursions, allowing, while seated, the control of exercise and its changes.

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This aid allows the control of movements between two segments. Thus, it is very easy to set different levels of flexion beyond which an error signal is activated (negative FB).

The patient is asked to achieve the maximum flexion, using a positive FB, until the indicator is activated.

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When lacking scapulothoracic stability during abduction of the upper limb the goniometer activates the vibratory feedback on the muscle that should be recruited .

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It helps control any compensation or lack of trunk stability during exercises which strengthen muscles, thus making exercises more accurate and effective.

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Flexor strengthening exercise of the upper limb. While in the previous exercise the inclinometer signaled compensations on the frontal plane now the patient should not stretch the trunk when using the elastic band.

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
The patient is involved in pointing tasks by head movements; the acoustic FB system signals insufficient pelvic stabilization. Use of the laser pointer is particularly useful due to the difficulties in quantifying head movements and asking the patient to perform specific tasks. For example, the patient is asked both to achieve a position with eyes closed and then to check whether the required task has been performed.

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Head rotation without involving the trunk is required. The movement of the antenna and further contact with the aid closes the circuit and activates the indicator. The device also helps quantify the movements and any asymmetry.

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For patients having difficulty straightening their trunks, a positive FB can be used until they return to their starting position or alternatively the laser can be used if they want a visible track of the trunk behaviour on the sagittal/horizontal plane.

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The inclinometer connected at the back of the thoracic indicator warns the patient that he/she has exceeded a predefined range of flexion. Positioning the vibrating pad on the dorsal muscle facilitates control of the proper backbone alignment.

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The patient is asked to identify stops placed laterally at different heights under the board while transferring the load from side to side without moving the upper part of the trunk. Any inclination (to the left in the photo) activates the FB.

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While the patient is seated at the oscillating board, the number of springs on either side is changed, thus influencing the load distribution.

The patient is asked to keep his/her trunk on axis while being assisted by the inclinometer which signals the movements exceeding the inclination on the frontal plane.

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The patient is asked to transfer his/her load, and to stop at the positions required by the physiotherapist. This exercise allows initial work involving load variation and distribution within the support base. As with the previous exercise, the inclinometer is placed on the trunk.

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With his/her hands placed on the board with spheres, the patient is asked to assume different positions on the sliding surface while keeping the pointer on the target at the front.

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
This exercise is particularly useful to define the change from sitting to standing position. In this case, it is important to provide the patient with a reference beyond which he/she is to recruit the muscles of the lower limbs to pull away his/her pelvis. The patient will be asked to continue moving to the next reference and back, while keeping his/her hands placed on the board.

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The patient is taught to bend forward at his/her hips.

The goniometer activates the indicator in case of excessive pelvic retroversion.

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In order to improve the knee joint function, it is better to apply the laser to the ankle as precise reference to define angle variations. The patient must follow the physiotherapist and identify different positions.

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The laser connected to the board magnifies the ankle movements also providing a visual feedback. When using the clasps the patient works on a single axis of rotation while the hemispheres enable the rotation exercises.

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The patient identifies or reaches different board positions on which he/she places a foot while performing tasks to flex/extend the knee. The inclinometer on the femur is activated when the limb rotates inward.

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The goniometer positioned at the front of the thigh extends during hip extension, thus activating a positive FB. By measuring the time the FB remains active it is possible to act to improve the recruitment.

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In static posture, with segments aligned, the fixed antenna can reproduce the movement of the barycentre. The oscillations of the body, consequently the barycentre of the patient, are graphically indicated by the tip of the antenna which keeps track on the coloured panel.

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
Sensory manipulation exercise. By reducing the information from the plantar surface and asking the patient to close his/her eyes, he/she is "forced" to find the barycentre position by the vestibular system. The antenna on the pelvis keeps tracks of it.

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This exercise is useful to control the barycentre. The patient is asked to flex/extend his/her trunk while on the board, which is stabilized by springs in front and blocked behind by a stop. Excessive transfer of his/her load on the forefoot activates the inclinometer.

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This exercise helps providing accurate and reliable information to patients suffering from balance disturbance and with difficulty to perceive centre-of-pressure variations.

The vibratory FB located on the posterior part of the leg is activated upon inclination of the leg beyond a pre-established range.

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
The patient allows load to be moved within the support base on the sagittal plane while using ankles as a fulcrum (inverted pendulum model).

The pelvic indicator allows precisely quantifying the extent of anteroposterior movement. If the light target of the laser pointer on the trunk moves excessively, it means the patient is not using the correct strategy for the ankles but also involves other joints.

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
Search for posterior Limits Of Stability (LOS). In order to avoid that patient's trunk leans forward, the inclinometer is placed on the back of his/her trunk. Backing off of antenna while the trunk is stabilized corresponds to a backward movement of the barycentre.

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
Exercise similar to the previous, with the difference that in this case the search is for anterior limit of the displacement of the barycentre. The inclinometer positioned at the front of the trunk signals undesired extension movements.

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The patient is asked to identify the number of washers located under the board while keeping visual contact with the pelvic indicator and ensuring his/her knee remains extended. The washers are placed on one of the contacts while the other is screwed into the board. During this exercise it is possible to obtain precise acoustic feedback at extreme levels of dorsal flexion, thus making the work of extending the posterior muscles effective.

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
The patient on the oscillating board is asked to lower his/her heel to identify the different heights under the board. The inclinometer connected to the tibia signals any possible backward movement of the tibia. In this manner, it is possible to work on tibiotarsal articulation in a selective way.

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The patient is asked to let him/herself be moved on the frontal plane and to perceive the final position. The use of a grid facilitates reading motor behaviour. It is possible to ask the patient to move by him/herself while keeping the antenna parallel to the reference lines. The same exercise is used also to guide movements on the sagittal plane.

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Thanks to numerous holes and using 2 magnetized pins and 2 clasps of different heights, it is possible to create oscillating boards suitable for the different exercises. By inserting a magnetized pin between the board and the clasp, not only is it possible to change the position of the fulcrum, but also the axis with respect to the foot. Inserting the pointer in the extremity of the clasp facilitates the visualization of this. At this point, the patient is asked to transfer his/her load and to continuously modify recruitment in relation to the position of the clasp. Springs facilitate control.

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This aid is similar to the previous, with the difference that in this case, foot rotation axis allows control movements on the frontal plane.

The patient is asked to identify the direction of the clasp while dynamically transferring his/her load.

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Once identified the position of the axis of rotation, the patient is asked to transfer the load to the inner or outer part of the foot in order to identify the number of springs that are positioned, from time to time, on the sides of the board.

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Unlike the previous exercise, the task here is to be able to identify the board inclination, which is stopped by cylinders of different heights.

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The physiotherapist guides the patient's leg on the sagittal plane and by the indicator he identifies references. The patient is asked to identify or achieve knee angular positions while keeping the antenna stable on the pelvis.

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Depending on the problem to solve, the patient may be asked to complete the knee extension or give a signal to stop so as to avoid hyperextension.

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The patient is asked to achieve angular positions on the horizontal plane using the antenna at the knee. The patients is asked to keep the light target of the laser pointer on the pelvis within the range indicated so as to allow performing correct hip rotation.

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The patient is asked to flex his/her knee, loading the limb and keeping a specific angle of rotation. The laser pointer, placed distally at the femur, projects information on the panel in relation to the extent of flex/extension along with a correct rotation.

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
The scale provides information on the load distribution. The leg ahead rests on the oscillating board; the laser shows the flexion at knee and the rotation of the limb. The difficulty of this exercise depends on the number of springs placed on the sides of the board

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
The patient controls and distributes load on his/her lower limbs, keeping his/her balance while standing on the board that is stabilized laterally with springs. The patient performs tasks to flex/extend the knee; the laser on the knee shows the flexion level, and, keeping the pointer to the right (see photo), the knee is guided so as not to work in internal rotation .

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Pelvis retroversion exercise. The patient is asked to activate the FB connected to the inclinometer while keeping the laser pointer on the trunk fixed on a reference point that is outlined on the panel.

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
The aim of this exercise is to control anterolateral pelvis movement when it is fully loaded (to the left in the photo). Curvilinear trajectories with an increasing radius of curvature are outlined on the panel, and generated from the same point. Each trajectory guides progress by the patient's hemiside in the flight phase and represents a different level of forward movement of the leg on the ground.

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Graphic representation of the feet on the reference plane helps the patient transfer the load with appropriate movement of the pelvis on the horizontal plane .

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The patient is helped to complete the dorsiflexion at ankle during the weight-bearing phase, asking him to wait for the acoustic signal from the goniometer before raising his/her heel from the ground.

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By positioning the inclinometer on the pelvis (as shown in the photo) we can control movements when raising the right hemipelvis or the Trendelenburg position to the left.

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The patient is asked to move the laser pointer's light reference from position A to position B by rotating his/her pelvis on the horizontal plane.

Attention is drawn to the behaviour by the hemipelvis of the limb being prepared for the flight phase. During movement, if the hemipelvis is raised, the circuit is closed and the acoustic indicator is activated.

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
Similar to the previous exercise but in this case, the patient's attention is also drawn to the control of the knee which has to achieve a specific flexion.

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The patient is asked to relax the limb, by flexing the right knee until the goniometer is activated in the flight phase.

The laser applied at the side of the pelvic indicator describes the behaviour of the pelvis on the frontal plane, so as to quantify the elevation of pelvis or trendelenburg.

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The aim of this exercise is the correct contact of the foot with the ground. To reach an optimal dorsiflexion at ankle, the FB is required to be activated in the final phase of flight. The positive FB can be also determined by the vibrating pad on the tibialis muscle.

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
The aim of this exercise is achieve correct timing of the plantar flexors during forward thrust. The goniometer applied to the back of the left hip is activated upon maximum flexion in the final flight phase, thus operating the vibrating pad that is placed on the calf of the contralateral limb.

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The use of the laser pointer on the pelvis during ambulation allows further magnifying rotations. During the phase when the heel is on the ground, the homolateral hemipelvis is at maximum rotation and the pointer projecting on the wall or floor has reached its maximum lateral movement, and then changes direction. A study of the resulting sinusoid highlights pathological behaviour.

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An effective method to affect the gait length and symmetry is to apply a goniometer to the ankles, using special supports. Ankle clips release automatically (safety mechanism) if the person widely exceeds the length.

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A micro-mouse connected to the pelvis via an antenna, allows reading pelvic movements and store them on the PC. If trunk movements are avoided (inclinometer for the sagittal plane), shifts coincide with the shift of the center of mass.

CD installazione ed esercizi
Installation CD and exercises

MAGILL R.A. Motor learning: concepts and applications, 1th edition. Wm. C. Brown Group; Dubuque, Iowa: pp.14; 1980.

WOOLLACOTT M.H, SHUMWAY-COOK A. Maturation of feedback control of posture and equilibrium. In: Motor development in children. Fedrizzi E, Avanzini G, Crenna P. John Libbey & Co.; London: pp.59-70; 1994.

SMITH KN. Biofeedback in strokes. The Australian Journal of Physiotherapy. 25:155-161; 1979.

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CARR J.H, SHEPHERD R.B. Movement science: foundation for physical therapy in rehabilitation. Aspen Publishers; Rockville, Maryland, 1987.

ADAMS J.A. Historical review and appraisal of research on the learning, retention, and transfer of human motot skills. Psychological Bulletin. 101:41-74; 1987.

WINSTEIN C.J. Knowledge of results and motor learning: implications for physical therapy. Physical Therapy. 71 (2):140-149; 1991.

SCHMIDT R.A, LEE T.D. Motor control and learning: a behavioral emphasis, 4th edition. Human Kinetics Publishers; Champaign, IL, 2005.

WULF G. Attention and motor skill learning. Human Kinetics; Champaign, IL, 2007.

ROELFSEMA P.R, VAN OOYEN A, WATANABE T. Perceptual learning rules based on reinforcers and attention. Trends in Cognitive Sciences. 14 (2):64-71; 2010.

JONSDOTTIR J, CRIPPA A, PARELLI R, RABUFFETTI M, CASIRAGHI A, CATTANEO D.

Positional biofeedback to improve gait in subjeet with stroke
Gait & Posture Volume 30, Supplement 1 October 2009

JONSDOTTIR J, CATTANEO D, RECALCATI M, REGOLA A, RABUFFETTI M, FERRARIN M, CASIRAGHI A, Task-Oriented Biofeedback to Improve Gait in Individuals With Chronic Stroke: Motor Learning Approach
Neurorehabil Neural Repair 2010 24: 478 originally published online 6 January 2010

JONSDOTTIR J, CRIPPA A, CATTANEO D, RECALCATI M, REGOLA A, RABUFFETTI M, FERRARIN M, CASIRAGHI A. Concepts of Motor Learning Applied to a Rehabilitation Protocol Using Biofeedback to Improve Gait in a Chronic Stroke Patient: An A-B System Study With Multiple Gait Analyses
Neurorehabil Neural Repair 2007 21: 190

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WINTER Biomechanics and motor control of human gait. Waterloo: University Press. 1987

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REHABILITATION and MEDICAL EQUIPMENT

CHINESPORT spa - Via Croazia, 2 - 33100 Udine - Italy
Tel. +39 0432 621 621 - Fax +39 0432 621 620

chinesport@chinesport.it - export@chinesport.it

www.chinesport.it

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© Chinesport spa 2012