T-one

EVO I

EVO II

**Manuale utente
User Manual
Mode d'emploi**

I.A.C.E.R. Srl

Via S. Pertini 24/A, 30030 Martellago (VE) Italy

Tel. +39 041 5401356 Fax +39 041 5402684

e-mail: iacer@iacer.it <http://www.iacer.ve.it>

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T-ONE

EVO I EVO II

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EVO I EVO II

**WARNING: T-ONE IS
A MEDICAL DEVICE.**

Consult your doctor before using T-ONE
if you suffer from health problems.

Read the user manual and the electrode position manual before using T-ONE. Read the contraindications and warnings carefully before using the medical programs.

IACER S.r.l. via Pertini 24/A, 30030 Martellago (Venezia)-Italy.
Tel. +39 041 5401356 Fax +39 041 540 2684
iacer@iacer.it www.itechmedicaldivision.com

1. Presentation	Page 31
2. T-ONE kit contents;	Page 31
3. Purpose Contraindications and side effects	Page 31
4. Warnings	Page 32
5. Legend of symbols marked on T-ONE and found in the manual	Page 33
6. How to use T-ONE Description of controls – Preliminary instructions – Operating instructions	Page 33
7. Battery charging	Page 36
8. T-ONE maintenance and storage	Page 36
9. Resolution of problems and questions	Page 37
10. Electrostimulation and T-ONE stimulation intensities	Page 37
11. Sport programs	Page 38
12. Beauty programs	Page 43
13. T.E.N.S. programs	Page 46
14. SETUP	Page 51
15. Self-adhesive electrodes	Page 52
16. Declaration of Conformity	Page 52
17. Technical specifications	Page 52
18. Assistance	Page 53

1. Presentation.

Thank you for choosing the T-ONE electronic stimulator; this excellent medical device has been developed using modern and innovative technology, designed by a leading Italian company in the design and manufacture of professional electronic stimulators and surface electromyographs.

T-ONE EVO I: 11 Programmi Sport, 9 Programmi Beauty, 5 Programmi TENS

T-ONE EVO II: 11 Programmi Sport, 9 Programmi Beauty, 16 Programmi TENS

2. T-ONE kit contents.

1 T-ONE electronic stimulator;
4 two connection cables, for the transmission of electric impulses;
2 pkgs. containing 4 pre-gelled self-adhesive 41x41 mm electrodes (1 for T-ONE EVO I);
2 pkgs. containing 4 pre-gelled self-adhesive 40x80 mm electrodes (1 for T-ONE EVO I);
1 pkg. containing 8 pre-gelled diam. 32mm self-adhesive electrodes (not available on T-ONE EVO I);
1 battery charger;
1 T-ONE user manual;
1 electrode position manual;
1 carry bag.

Accessories:

- Pre-gelled electrodes: 41x41mm • 48X48mm • 40x80 mm • 50x90mm • diam. 32mm;
- Reusable elastic belts for limbs and calves;
- Electrotherapy gel

3. Purpose.

The T-ONE electronic stimulator is a medical device specifically intended for adults with adequate knowledge of the device acquired from reading the manual for domestic use.

T-ONE is used to apply electrical microimpulses, capable of producing energy, to the human body. This energy, modulated in accordance with parameters typical to each impulse, can be used for various purposes, from pain reduction to winding-down exercises, from muscle strengthening to drainage, from isotonic exercises to treating haematomas, to the treatment of imperfections using the beauty programs.

The manufacturer recommend to consult your doctor before using TENS programs.

Contraindications.

This device must not be used by people with a pace-maker, epileptics, pregnant women, people with a high temperature, anxious people or people suffering from heart disease, phlebitis, thrombophlebitis or serious illnesses.

In the event of injury, muscle stress or any other health problem consult your doctor before using the device and only use it under medical supervision.

Side effects

No significant side effects are known of. Some particularly sensitive people have reported skin redness in the area where the electrodes were positioned: the redness usually disappears a few minutes after the end of the treatment. Should the redness persist please consult a doctor. In rare cases, stimulation carried out in the evening can cause some people to experience difficulty in falling asleep. If this occurs, suspend the treatment and consult a doctor.

4. Warnings.

T-ONE has been designed and manufactured to be operated exclusively with the internal, rechargeable battery.

- T-ONE doesn't produce and doesn't receive electromagnetic interferences from other equipments. However, it's opportune to use the instrument keeping it at least at 3 meters from televisions, monitors, cellular telephones or any other electronic equipment.
- The device must not be used in the presence of patient monitoring equipment.
- Do not use T-ONE with electrosurgical or shortwave or microwave therapy equipment .
- Use of the device is prohibited to persons known to be of unsound mind, to persons suffering from sensibility disorders and to permanently or temporarily disabled persons unless assisted by qualified personnel (e.g. a doctor or therapist).
- Use of the device is prohibited in the presence of signs of deterioration of the device itself.

Should any foreign materials penetrate the device contact the retailer or manufacturer immediately.

If dropped down, check that the housing is not cracked or damaged in any way; if so, contact the retailer or manufacturer.

Should you notice any changes in the device's performance during treatment, interrupt the treatment immediately and consult the retailer or manufacturer (patients being treated in a centre must also be informed of the event).

- Use of the device is prohibited in combination with other medical devices.
- Use of the device is prohibited close to flammable substances or in environments with high concentrations of oxygen.
- Consult a doctor before using T-ONE with metallic osteosynthesis devices.

Muscle heaviness and cramps may be experienced in the event of prolonged muscle stimulation. Should this occur we recommend to suspend treatment with T-ONE for a few days and then to continue with muscle recovery and capillarization programs.

Make sure that the plug and battery charger are intact. Should one of these parts show any sign of damage, suspend use of the battery charger immediately and contact the retailer or manufacturer.

- It is forbidden to position the electrodes in such a way that the current crosses the heart area (e.g. a black electrode on the chest and a red electrode on the shoulder blade); however electrodes can be positioned along the muscular fascia of the heart area, as used for pectoral strengthening.
- Use of the device is prohibited with electrodes positioned on or close to injuries or cuts.
- The electrodes must not be positioned on the carotid sinuses (carotid) or genitals.
- The electrodes must not be positioned close to the eyes; make sure that the current delivered does not cross the eyeball (one electrode diametrically opposite to the other in relation to the eye); keep a distance of at least 3 cm. from the eyeball.
- When using the electrodes, follow the instructions given in this manual and on the package of the electrodes. Use single-patient electrodes, supplied exclusively by the manufacturer, and take care to avoid the exchange of electrodes between different users. T-ONE has been tested and guaranteed for use with the electrodes supplied.
- Insufficiently sized electrode sections can cause skin reactions or burns. Only use the round electrodes supplied by the manufacturer for facial applications.
- Do not use electrodes when damaged, even if they stick to the skin well.
- Electrodes must not be used when they no longer stick to the skin. Repeated use of the same electrodes can compromise the safety of the stimulation, in fact it can cause skin redness that can last for many hours after stimulation.
- Warning! During treatment with round and 35 x 45 mm. electrodes, current density may be in excess of 2 mA/cm² for each electrode. In these cases take special care regarding any skin redness.

Product expected lifetime: 5 years

5. Legend of symbols marked on T-ONE and found in the manual.

Produced in compliance with the 93/42/EEC directive concerning medical devices

Internally powered device with type BF applied parts

Class II equipment

Electrostatic discharge sensibility (ESD)

Product subject to WEEE regulations concerning separate waste collection

Read the instructions for use carefully before using the device

ON/OFF button.

Select program and muscles to stimulate.

Select Sport, Beauty, Tens programs, male and female.

Start/Stop program, confirm selection.

Increase intensity.

Decrease intensity.

6. How to use T-ONE.

Preliminary instructions.

Before using T-ONE clean the skin of the area to be treated; with the cable disconnected from T-ONE, connect the electrostimulation cable jacks to the self-adhesive electrodes; position the self-adhesive electrodes on the skin (see photos of electrode positions); connect the impulse transmission cables to the relative jacks, then turn T-ONE on.

Warning: make sure that T-ONE is switched off before disconnecting the electrodes at the end of the treatment.

Operating instructions (N.B. read the whole manual before using the equipment)

1. Turn T-ONE on using the  button
2. Select the program group (Sport, Beauty, Tens) using the button  
3. Select the program using the   buttons


Sport programs display




Beauty programs display


4. Select male or female with   buttons .

5. Select the muscle group with   buttons

6. press  to start program; the display shows connect cables and electrodes and increase intensity of the used channels; to start intensity press  for each used channel; T-ONE will begin the session.
7. press the  button for each use channel to increase current intensity until the personal tolerance level is reached.

The picture above is only used as demonstration and it not represent exactly the specifications of program S10 Endurance

A bar indicating the start and end of a contraction will appear on the display whilst the electronic stimulator is working. The periods of muscular contraction (bold bar) and the periods of recovery (striped bar) are graphically visualised under the bar with a cursor. This enables the user to know exactly when a contractions starts.

8. At the end of the first phase, the T-ONE cancels the previously selected intensity, and warns the user with an intermittent signal; to continue with the program, press the button for the relevant channels again until you reach the desired level of stimulation. An audio signal will sound three times when the program has finished.
9. At the end of the program switch off the T-ONE and remove the cables. Attach the electrodes to the transparent films and replace them in their original packaging.

Option to simultaneously increase the intensity on all channels:

Select the desired program:

Start program

Increase each channel to at least intensity 1:

CH1
 CH2
 CH3
 CH4
 Blue Yellow Red White

To simultaneously increase the intensity of all channels use the following buttons:

Skip phase: press the  button during the working program to skip to the next phase.

Pause/stop program: during the working program, press the  button once to momentarily pause the program. Press the  button again to restart the program.
Switch off the machine to stop the program.

7. Battery charging.

The display shows when the battery is low of power. Insert the jack of the battery charger into its socket and plug the battery charger into the main power socket.
Recharge for 8 hours for a complete charge.

Do not immerge the battery in water or other liquids and do not expose it to heat sources.
Do not dispose of dead or defective batteries with domestic waste; dispose of in an authorised waste collection bin.

Before storing T-ONE for long periods of time, we recommend that it is fully charged and stored in a dry place, checking that it is switched off.

Do not use the battery charger if:

- The plug is damaged or parts of it are broken.
- It has been exposed to rain or any other type of liquid.
- The components have been damaged by a fall.

Do not immerge the battery charger in water. Use a dry cloth to clean the battery charger. Do not open the battery charger: it does not contain repairable parts.

Warning: disconnect electrodes during recharchig.

8. T-ONE maintenance and storage.

T-ONE must be switched off and the cables removed from the connectors at the end of every session. T-ONE should be kept in the carry bag supplied, together with its accessories, and carefully stored on a secure surface and protected from situations listed in the *Warnings* paragraph.

T-ONE should be stored in the following conditions:

Temperature: from -20°C to +40°C
Maximum relative humidity: 75%
Atmospheric pressure: from 700 hPa to 1060 hPa

Cleaning: only clean the device using a dry soft cloth.

We recommend a functional test carried out by the manufacturer every 24 months.

The manufacturer retains that the T-ONE device cannot be repaired by personnel external to its own company. Any work of this nature carried out by personnel not authorised by the manufacturer will be classed as tampering with the device, freeing the manufacturer from responsibility for the warranty and hazards that the operator or user may be exposed to.

N.B. disconnect the cables before placing the electronic stimulator in the bag. If not, the cables will bend excessively close to the connectors. This may damage the cables.

9. Resolution of problems and questions.

Any type of work on T-ONE must be carried out exclusively by the manufacturer or by an authorised dealer. In any event, any presumed malfunction of T-ONE must be verified before sending the device to the manufacturer.

Here below are some typical situations:

T-ONE cannot be turned on:

- Recharge the batteries and try to switch it on again; if it still does not switch on, check that the power supply socket is working.
- Check that the ON/OFF button has been pressed correctly (held down for at least one second).

T-ONE does not transmit electric impulses:

- Check that the cable jacks have been inserted in the electrodes and that the plastic protection has been removed from the electrode.
- Check that the cables have been connected correctly (connector well inserted in the device).

T-ONE switches off during operation:

- Recharge the batteries for 8 hours and start the treatment again. If the problem persists contact the manufacturer.

T-ONE is switched on but nothing appears on the display:

- Open the battery compartment, lift up the battery, disconnect the connection and remove it. Inside battery compartment you'll find a little hole with RESET writing above: insert a little tip and gently press.

Disposal

To protect the environment, the device and its accessories must be disposed of in accordance with the law in areas equipped for that purpose or as special waste.

Product subject to WEEE regulations concerning separate waste collection: when disposing this product, please use the designated areas for disposing electronic waste or contact the manufacturer.

10. Electrostimulation and stimulation intensity.

Electrostimulation consists of the transmission of electric microimpulses to the human body. Fields of application of electrostimulation are: pain therapy, recovery of muscle trophism after injury or a surgical operation, athletic preparation and beauty treatments. Specific electric impulses are used for every one of these applications.

The stimulation intensity is displayed on the T-ONE screen for each channel in a scale progressing from 0 to 120.


Once the initial intensity has been set, T-ONE will proceed until it reaches the end of the phase; when more than one phase is present, the intensity has to be set again at the beginning of new phase.

The types of impulses can be classified as follows:

1. **Tens impulse:** for tens programs, the intensity should be adjusted to a level between the thresholds of perception and pain. The maximum limit is reached when the muscles surrounding the area treated begin to contract. We recommend to stay below that limit.

2. **Lipolysis/drainage impulse:** the “pump” effect is produced by sequential tonic contractions. The intensity must be enough to produce these contractions: the greater the contraction is, the greater the pump effect will be. But beware: there is no benefit to be gained from enduring intensity high enough to cause pain. The first electrostimulation sessions should be carried out at a low intensity to allow the organism to get used to new sensations. That way intensity can be increased gradually and will not be traumatic.
3. **Warming-up impulse:** stimulation intensity must be increased gradually to improve the metabolism of the muscle concerned progressively. A bit like with a car: the engine needs warming-up before it is taken to full speed.
4. **Toning, training, atrophy contraction impulse:** the muscle treated must visibly contract during a training impulse. The fact that the muscle tends to stiffen and increase in volume will be visible. Intensity should be increased gradually (in the first contraction) to enable you to identify the right level of stimulation comfort. Intensity can be increased up to the personal tolerance threshold during the second training contraction; this operation is then repeated during each contraction until the workload reaches the level of intensity recommended in the description of the single programs. We recommend that you record the level of intensity reached in order to improve the level of stimulation and consequently your performance.
5. **Massage, winding down, active recovery impulse:** intensity must be adjusted gradually to massage the muscle treated. The level of intensity should be sufficient to obtain a comfortable massage. There is no need to endure high intensity levels in this case as it is meant to be a massage, meaning that intensity can be increased gradually without excesses.
6. **Capillarization impulse:** increase the intensity gradually to produce constant, visible stimulation of the area treated; a medium stimulation threshold is recommended, always below the pain threshold.

11. Sports programs.

 **IMPORTANT. Stimulation intensity:** the sports programs (Resistance, Strength, Agonist/Antagonist) are divided into warm up phases, work phases, and recovery phases. During the warm up phase (phase 1), the stimulation intensity should be adjusted to produce a slight tickling sensation in the treated muscle, this should be enough to warm up the muscle without tiring it (18÷30 mA). During the work phases (the central phases of the programs), contraction impulses are alternated with active recovery impulses. The user needs to select the two intensities separately: for the contraction, the user should see the following on the display with a bold bar; 20÷30 mA for relatively unfit people, 30÷50 mA for fit people, and more than 50 mA for very fit people.

We recommend adjusting the intensity of the recovery phase between one contraction and another by approximately 15÷20% less than the intensity of the contraction. During the recovery phase (the last phase), the intensity should be adjusted to produce a good massaging effect of the stimulated area without producing any pain (18÷30 mA).

S01 • Warming up (all muscle groups).

Program suitable for use before training sessions or competitions, very useful for sports involving maximum effort right from the start. Program duration: 28 minutes. Position of electrodes from photo 1 to photo 20. Recommended stimulation **intensity:** medium (20-30mA); the muscle must work without strain.

S01 Warming up		Total duration: 28 minutes
	Biceps (photo 02/15), Triceps (photo 03/16), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).	
	Abdominals (photo 01/20), Pectoral/Breast (photo 07/17), Trapezius (photo 08), Large dorsal (photo 09), Gluteus (photo 19)	
	Quadriceps/Thighs (photo 11/18), Biceps femoris (photo 12), Calves (photo 13), Anterior Tibial (photo 14).	
Phase 1	Phase 2	Phase 3
3 min. Warm up with low frequency	5 min. Warm up with medium frequency	20 min. 20 sec. recovery 10 sec. contraction

S02 • Recovery.

Can be used for all sports, after competitions or the most demanding training sessions, in particular after long and intense exertion. To be used immediately after exertion. Helps drainage and winding down, improving muscle oxygenation and helping to discharge synthetic substances produced during exertion. Program duration: 15 minutes/3 phases. Stimulation **intensity**: medium-low, increased during the last 5 minutes.

Position of electrodes: **see photos 01 to 20 in the manual of positions in relation to the area that you wish to stimulate.**

S03 • Capillarization.

The capillarization program significantly increases arterial flow in the area treated. The capillarization program is very useful for recovering after intense aerobic work (toning up training) and improves local microcirculation. Program duration: 30 minutes. Recommended stimulation **intensity**: medium (20-30mA). Position of electrodes: see photos 01 to 20 in the manual of positions.

S04 • Competitive/Counteractive (Agonist / Antagonist)

The electronic stimulator produces contractions alternated between 4 channels: during the first 3 minutes of warm-up the 4 channels work simultaneously, during the second work phase (20 minutes) muscle contractions are alternated between Channel 1/2 (agonist muscles) and Channel 3/4 (antagonist muscles). The program is designed to restore muscle tone to the quadriceps and its antagonist the leg biceps, or the biceps brachii and the triceps. The work aims at developing strength. With this program, muscle relaxation is obtained by simultaneous stimulation from 4 channels during the last 10 minutes (phase 3).

Stimulation **intensity** during the contraction: enough to produce good muscle contraction + voluntary contraction to reduce the sense of discomfort and reach higher intensities. Intensity must be increased gradually treatment by treatment, without overstraining the muscles. Suspend training for a few days in the event of fatigue and proceed with the "S 02" program.

During the second work phase, the intensity level needs to be selected three times: once for the four channels of active recovery, once for contractions with channels 1 and 2 (agonist muscles) and once for contractions with channels 3 and 4 (antagonist muscles).

S04 Agonist / Antagonist		Total duration: 33 minutes
	Biceps brachii (CH1+CH2 - photo 02) / Triceps (CH3+CH4 - photo 03), Hand extensors (CH1+CH2 - photo 04) / Hand flexors (CH3+CH4 - photo 05).	
	Quadriceps (CH1+CH2 - photo 11) / Biceps femoris (CH3+CH4 - photo 12).	
Phase 1	Phase 2	Phase 3

3 min. warm up	20 min. alternated contractions on 4 channels 6 sec. Recovery CH1/2/3/4 7sec. contraction CH1/CH2 7sec. contraction CH3/CH4	10 min. recovery
----------------	--	------------------

S05 • Long effort (resistant strength).

This program is designed to increase resistance to physical stress, or rather withstand intense exertion for a longer amount of time in muscle regions subjected to stimulation. Indicated for sporting disciplines involving long and intense periods of exertion.

Stimulation **intensity** during the contraction: start with a low level of intensity, increasing it gradually. For practised athletes the intensity used should be enough to produce visible muscle contractions.

In the event of muscle ache after stimulation, use the S02 program (recovery).

S05 Long effort		Total duration: 30 minutes	
	Biceps (photo 02/15), Triceps (photo 03/16), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/20), Pectoral/Breast (photo 07/17), Trapezius (photo 08), Large dorsal (photo 09), Gluteus (photo 19)		
	Quadriceps/Thighs (photo 11/18), Biceps femoris (photo12), Calves (photo13), Anterior Tibial (photo14).		
Phase 1	Phase 2	Phase 3	Phase 4
5 min. warm up	10 min. 20 sec. recovery 10 sec. contraction	10 min. 20 sec. recovery 10 sec. contraction	5 min. recovery

S06 • Basic strength.

The Basic strength program is used in sport to develop basic strength, which for definition is the maximum tension that a muscle can exert against constant resistance. The contractions alternate with periods of active recovery during the work phase, allowing the muscle to be trained without subjecting it to stress and improving oxygenation of the same muscle.

The following basic procedure will enable you to obtain the first results: two sessions per week (for each muscle region) for the first three weeks at medium/low intensity (20-30mA), three sessions per week for the next three weeks at high **intensity** (30-50mA).

Intensity must be increased gradually treatment by treatment, without overstraining the muscles. Suspend training for a few days in the event of fatigue and proceed with the "S 02" program.

S06 Basic strenght		Total duration: 30 minutes	
	Biceps (photo 02/15), Triceps (photo 03/16), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/20), Pectoral/Breast (photo 07/17), Trapezius (photo 08), Large dorsal (photo 09), Gluteus (photo 19)		
	Quadriceps/Thighs (photo 11/18), Biceps femoris (photo12), Calves (photo13), Anterior Tibial (photo14).		
Phase 1	Phase 2	Phase 3	Phase 4
5 min. warm up	10 min. 13 sec. recovery 7 sec. contraction	10 min. 14 sec. recovery 6 sec. contraction	5 min. recovery

S07 • Speed (Fast strength).

This program is designed to increase speed in fast athletes and develop it in athletes lacking the quality.

The exercise assumes a fast pace and the contraction is short, as is the recovery. It is usually best to complete a three-week basic strength cycle of increasing intensity before using this program. Then continue with three weeks of fast strength three times a week at high **intensity** (30-50mA) during the contraction.

S07 Fast strength		Total duration: 26 minutes	
	Biceps (photo 02/15), Triceps (photo 03/16), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/20), Pectoral/Breast (photo 07/17), Trapezius (photo 08), Large dorsal (photo 09), Gluteus (photo 19)		
	Quadriceps/Thighs (photo 11/18), Biceps femoris (photo 12), Calves (photo 13), Anterior Tibial (photo 14).		
Phase 1	Phase 2	Phase 3	Phase 4
3 min. warm up	10 min. 14 sec. recovery 6 sec. contraction	8 min. 10 sec. recovery 5 sec. contraction	5 min. recovery

S08 • Explosiveness (explosive strength)

Explosive strength programs increase the explosive power and speed of the muscle mass, with extremely short, strengthening contractions and very long active recovery times to allow the muscle to regain strength. Complete a three-week basic strength cycle (S06) before using this program. Then continue with three weeks of explosive strength twice a week. During the contraction, the **intensity** must be the highest that can be endured in order to obtain maximum muscle exertion whilst involving the greatest number of fibres (over 35mA).

S08 Explosive strength		Total duration: 28 minutes	
	Biceps (photo 02/15), Triceps (photo 03/16), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/20), Pectoral/Breast (photo 07/17), Trapezius (photo 08), Large dorsal (photo 09), Gluteus (photo 19)		
	Quadriceps/Thighs (photo 11/18), Biceps femoris (photo 12), Calves (photo 13), Anterior Tibial (photo 14).		
Phase 1	Phase 2	Phase 3	Phase 4
3 min. warm up	10 min. 24 sec. recovery 6 sec. contraction	10 min. 24 sec. recovery 5 sec. contraction	5 min. recovery

S09 • Endurance.

The Endurance program is used in sports to increase muscle resistance, acting mainly on slow-twitch fibres.

Program indicated for endurance sports: marathon runners, cross-country skiers, ironman, etc. Program duration: 34 minutes. Stimulation **intensity** during the contraction: if not particularly fit, start with a low intensity then increase gradually. For practised athletes the intensity used should be enough to produce visible muscle contractions. In the event of muscle ache after stimulation, use the Fitness 19 program (muscle relaxant).

S09 Endurance		Durata totale programma: 50 minuti	
	Biceps (photo 02/15), Triceps (photo 03/16), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/20), Pectoral/Breast (photo 07/17), Trapezius (photo 08), Large dorsal (photo 09), Gluteus (photo 19)		
	Quadriceps/Thighs (photo 11/18), Biceps femoris (photo 12), Calves (photo 13), Anterior Tibial (photo 14).		
Phase 1	Phase 2	Phase 3	Phase 4
5 min. warm up	20 min. 20 sec. recovery 10 sec. contraction	20 min. 20 sec. recovery 10 sec. contraction	5 min. recovery

S10 • Tonic contractions (sequential).

This program increases microcirculation within and around the muscle fibres treated creating rhythmic contractions, fostering better drainage and toning. It can also be applied to older people to improve blood and lymphatic circulation in the lower limbs (e.g. applying CH1 to the right calf, CH2 to the right thigh, CH3 to the left calf, CH4 to the left thigh).

This program can be carried out using self-adhesive electrodes. Stimulation **intensity** must be sufficient to produce good muscle contractions during the treatment but not enough to cause any soreness. It mainly works on slow-twitch fibres.

S10 Tonic contractions		Total duration: 25 minutes	
Muscles: all the muscles you want to stimulate in a sequence			
Phase 1	Phase 2	Phase 3	Phase 4
3 min. warm up	10 min. 6 sec. Sequential contractions on the 4 channels	10 min. 6 sec. Sequential contractions on the 4 channels	2 min. recovery

S11 • Phasic contractions (sequential).

This program produces rhythmic contractions with a stimulation frequency typical of fast-twitch fibres. It is suitable for increasing muscle strength sequentially.

The program produces sequential phasic contractions on 4 channels. Stimulation **intensity** must be sufficient to produce good muscle contractions during the treatment but not enough to cause any soreness.

Unlike the previous program, this one uses a higher stimulation frequency during the contraction phase and therefore works mainly on fast-twitch fibres.

S11 Phasic contractions		Total duration: 25 minutes	
Muscles: all the muscles you want to stimulate in a sequence			
Phase 1	Phase 2	Phase 3	Phase 4
3 min. warm up	10 min. 6 sec. Sequential contractions on the 4 channels	10 min. 6 sec. Sequential contractions on the 4 channels	2 min. recovery

12. Beauty programs.

E01 • Lipolysis.

This specific drainage program increases microcirculation within and around the muscle fibres treated and creates rhythmic contractions, facilitating the discharge of algogenic substances and promoting lymphatic activity. It can also be applied to older people to improve blood and lymphatic circulation.

The program produces sequential tonic contractions (before CH1/2, after CH3/4), reproducing the typical effect of electronic lymphatic drainage.

There are no real limits of application for this program, which can be practiced until the desired result has been achieved.

Stimulation **intensity** must be sufficient to produce good muscle contractions during the treatment but not enough to cause any soreness. Duration 30 minutes.

The first results can usually be seen after 3/4 weeks practicing 4/5 sessions a week.

The programme produces sequential tonic contractions in the central work phase (first CH1/CH2 then CH3/CH4): this is known as electronic lymphatic drainage. The electrodes should be applied as follows: CH1/CH2 on the extremities of the limbs (e.g. calf or forearm) and CH3/CH4 on the upper part (e.g. thigh or biceps brachii).

Stimulation intensity: this should be enough to produce a good level of stimulation (without pain) in the treated areas (15+30 mA).

E01 Lipolysis		Total duration: 30 minutes
Upper muscles:	Hand extensors CH1 (photo 04 with 2 electrodes) / Hand flexors CH2 (photo 05 with 2 electrodes); biceps CH3 (photo 02 with 2 electrodes) / Triceps CH4 (photo 03 with 2 electrodes)	
Lower muscles:	Calves CH1 (photo 13 con 2 elettrodi) / Anterior tibial CH2 (photo 14 with 2 electrodes); Quadriceps CH3 (photo 11 with 2 electrodes) / Bicipes femoris CH4 (photo 12 with 2 electrodes)	
Phase 1	Phase 2	Phase 3
5 min. warm up	20 min. alternated contractions: 6 sec. recovery on 4 channels 7sec. contraction CH1/CH2 7sec. contraction CH3/CH4	5 min. recovery

E02• Refirming.

Indicated for firming up muscles; working mainly on slow-twitch fibres. Suitable for those who have never done any physical activity or have been inactive for a long period of time. Method of use:

1. identify the muscle to be treated. To obtain good results we suggest to treat just a few muscles at a time and complete the process described below;
2. position the electrodes as shown in the photo;
3. increase the intensity until the impulse can be felt (use a low intensity for the first session to help you to understand how the machine works);
4. during the program and over the next few days the intensity should be increased gradually so that muscle contractions are not painful;
5. during contraction generated by the unit, contract the muscle voluntarily;
6. a cycle of 15/20 applications must be completed before the first results can be seen; one application for each muscle every two days with a day of rest in between.

It is possible to work on pairs of muscles, for example thighs and abdominal muscles, treating thighs one day and abdominal the next day. Working on too many muscles at the same time is not recommended.

E02 Refirming		Total duration: 30 minutes	
	Biceps (photo 02/15), Triceps (photo 03/16), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/20), Pectoral/Breast (photo 07/17), Trapezius (photo 08), Large dorsal (photo 09), Gluteus (photo 19)		
	Quadriceps/Thighs (photo 11/18), Biceps femoris (photo 12), Calves (photo 13), Anterior Tibial (photo 14).		
Phase 1	Phase 2	Phase 3	Phase 4
5 min. warm up	10 min. 13 sec. recovery 7 sec. contraction	10 min. 13 sec. recovery 7 sec. contraction	5 min. recovery

E03 Tonification.

Indicated for toning up muscles, working mainly on fast-twitch fibres. Suitable for people who already practice moderate physical activity. Method of use:

1. identify the muscle to be treated. To obtain good results we suggest to treat just a few muscles at a time and complete the process described below;
2. position the electrodes as shown in the photo;
3. increase the intensity until the impulse can be felt (use a low intensity for the first session to help you to understand how the machine works);
4. during the program and over the next few days the intensity should be increased gradually so that muscle contractions are not painful;
5. during contraction generated by the unit, contract the muscle voluntarily;
6. a cycle of 15/20 applications must be completed before the first results can be seen; one application for each muscle every two days with a day of rest in between.

It is possible to work on pairs of muscles, for example thighs and abdominal muscles, treating thighs one day and abdominal the next day. Working on too many muscles at the same time is not recommended.

E03 Tonification		Total duration: 34 minutes	
	Biceps (photo 02/15), Triceps (photo 03/16), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/20), Pectoral/Breast (photo 07/17), Trapezius (photo 08), Large dorsal (photo 09), Gluteus (photo 19)		
	Quadriceps/Thighs (photo 11/18), Biceps femoris (photo 12), Calves (photo 13), Anterior Tibial (photo 14).		
Phase 1	Phase 2	Phase 3	Phase 4
5 min. warm up	12 min. 13 sec. recovery 7 sec. contraction	12 min. 13 sec. recovery 7 sec. contraction	5 min. recovery

E04 • Definition.

Indicated for defining muscles, working mainly on explosive fibres. Suitable for people who already practise good physical activity and wish to define their muscles in greater detail. Method of use:

1. identify the muscle to be treated. To obtain good results we suggest to treat just a few muscles at a time and complete the process described below;

2. position the electrodes as shown in the photo;
3. increase the intensity until the impulse can be felt (use a low intensity for the first session to help you to understand how the machine works);
4. during the program and over the next few days the intensity should be increased gradually so that muscle contractions are not painful;
5. during contraction generated by the unit, contract the muscle voluntarily;
6. a cycle of 15/20 applications must be completed before the first results can be seen; one application for each muscle every two days with a day of rest in between.

It is possible to work on pairs of muscles, for example thighs and abdominal muscles, treating thighs one day and abdominal the next day. Working on too many muscles at the same time is not recommended. Stimulation intensity: 20÷35mA.

E04 Definition		Total duration: 25 minutes	
	Biceps (photo 02/15), Triceps (photo 03/16), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).		
	Abdominals (photo 01/20), Pectoral/Breast (photo 07/17), Trapezius (photo 08), Large dorsal (photo 09), Gluteus (photo 19)		
	Quadriceps/Thighs (photo 11/18), Biceps femoris (photo 12), Calves (photo 13), Anterior Tibial (photo 14).		
Phase 1	Phase 2	Phase 3	Phase 4
5 min. warm up	10 min. 13 sec. recovery 7 sec. contraction	5 min. 14 sec. recovery 6 sec. contraction	5 min. recovery

E05 • Limbs modelling.

Due to a combination of capillarizing and toning impulses, this program helps mobilise fat in areas where it tends to accumulate. The electrodes should be positioned to form a square around the limbs area and can be applied daily using a medium intensity.

Program duration: 30 minutes / one phase. Recommended stimulation **intensity**: medium.

Position of electrodes: see photos 01 to 20 and photo 22 and 23 in the manual of positions.

E06 • Microlifting.

The following program, with a duration of 15 minutes, is used to tone facial muscles using a special impulse to improve both the appearance and the dynamism of facial muscles.

The position of the electrodes is shown in the manual of electrode positions (photo 24).

N.B. a minimum distance of 3 cm. must be kept between the electrode and the eyeball.

⚠ IMPORTANT: take care when adjusting the intensity as facial muscles because they are particularly sensitive; intensity should be increased gradually, starting with a very low level of stimulation (just above perception) and increasing with care until a good level of stimulation is reached, represented by good muscle activation.

⚠ IMPORTANT: it is not necessary to reach levels of intensity capable of causing discomfort! The equation “more pain = more gain” is completely misleading and counterproductive. Great and significant results are obtained through consistency and patience.

E07 • Relaxing massage.

Can be used for all sports, after competitions or the most demanding training sessions, in particular after long and intense exertion. Program similar to the previous one: however it uses lower frequencies with a greater capacity for vascularization. To be used immediately after exertion. Helps drainage and capillarization, improving muscle oxygenation and helping to discharge synthetic substances produced during exertion. Program duration: 15 minutes. Stimulation **intensity**: medium-low, increased during the last 10 minutes.

Position of electrodes: see photos 01 to 20 in the manual of positions.

E08 • General recovery.

Can be used for all sports, after competitions or the most demanding training sessions, in particular after long and intense exertion. To be used immediately after exertion. Helps drainage and capillarization, improving muscle oxygenation and helping to discharge synthetic substances produced during exertion. Program duration: 15 minutes. Stimulation **intensity**: medium-low, increased during the last 10 minutes.

Position of electrodes: see photos 01 to 20 in the manual of positions.

E09 • Drain up.

These specific drainage programs increase microcirculation within and around the treated muscle fibres and create rhythmic contractions, facilitating the discharge of algogenic substances and promoting lymphatic activity. They can also be applied to older people to improve blood and lymphatic circulation.

The programs produce sequential tonic contractions in phases 2 and 3 that reproduce the typical effect of electronic lymphatic drainage. However in this case, unlike the lipolysis program, the contractions are sequential: first CH1 then CH2, CH3, and CH4.

There are no real limits of application for these programs, which can be practiced until the desired result has been achieved.

Stimulation intensity must be sufficient to produce good muscle contractions during the treatment but not enough to cause any soreness (20÷30 mA).

The first results can usually be seen after 3/4 weeks of practicing 4/5 sessions a week

E09 Drain up			Total duration: 25 minutes
Muscles: All the muscles you want to stimulate in a sequence. Connect a channel to each muscle that you want to stimulate.			
Phase 1	Phase 2	Phase 3	Phase 4
3 min. warm up	10 min. Alternated contraction on 4 channel, 6 sec each	10 min. Alternated contraction on 4 channel, 6 sec each	2 min. recovery

13. T.E.N.S. / rehabilitation programs.

WARNING!

It must be remembered that an electronic stimulator is a very effective analgesic instrument and that pain can indicate various types of medical conditions!

The Tens1 to Tens20 programs described in this paragraph are analgesic. You are advised to read the manual carefully before using T-ONE.

**REMEMBER THAT PAIN IS A SIGN:
consult a doctor to identify the medical condition before using T-ONE!**

TENS, an acronym standing for “transcutaneous electrical nerve stimulation”, is a therapeutic technique mainly used for analgesic purposes to counter the effects (usually pain) of a wide variety of medical conditions. For this purpose it finds application in treating everyday ailments troubling mankind: neck pain, arthrosis, myalgia, neuritis, back pain, periarthritis, heaviness in legs, muscle weakness, just to mention a few.


On an academic level, TENS can be divided into various categories according to the mechanism used to reduce the pain. The main types are: conventional TENS (or fast

analgesic), electro acupuncture training TENS (or delayed analgesic), TENS at maximum values with antidromic action and consequently an immediate local anaesthetic effect.

The rehabilitative action of TENS is represented by its power to reduce pain thereby restoring physiological conditions; this allows the patient to regain normal motor function. Consider a patient suffering from irritating peri-arthritis; they either resort to using analgesics or learn to live with the pain that often makes even the simplest movements impossible. Immobility reduces metabolic activity making it impossible to eliminate algogenic substances. So a vicious circle begins. In addition to relieving pain, TENS causes induced muscle stimulation increasing metabolic activity and blood flow and improving tissue oxygenation with an intake of nutritional substances. So the positive effect can be amplified by combining TENS with muscle stimulation of the area concerned.

Position of electrodes and intensity levels.

The electrodes should be positioned to form a square surrounding the painful area using Channel 1 and Channel 2 as shown in illustration 1 (red up/black down are not important for therapeutic purposes, follow instructions on electrodes position manual). The intensity should be adjusted to a level between the thresholds of perception and thresholds of pain: the maximum intensity level is the moment in which the muscles on the area treated begin to contract; over this limit, stimulation does not become more effective, just more irritating.

A01 • TENS rapid (Medical program)

Program also called **conventional tens**, used for analgesic purposes, blocking pain at the dorsal horn of the spinal cord, in accordance with the “Gate Control Theory” by Melzack and Wall. Pain impulses leave part of the body (for example the hand) and run along the nerve tracts (through small-diameter nerve fibres) until they reach the central nervous system where the impulses are interpreted as pain. Conventional tens activates large-diameter nerve fibres, blocking the path of small-diameter nerve fibres at the spinal cord. So action is mainly taken against the symptom: to simplify it further, the wire transmitting pain information is obstructed. Treatment duration should be no less than 30/40 minutes. **Conventional tens is a current that can be used to treat general everyday pain.** The average number required to benefit from the treatment is 10/12 once a day (there are no contraindications for up to double this amount).

The program has a duration of 30 minutes in a single phase. The program can be repeated at the end of the session for particularly persistent pain. The nature of the impulse means that the patient may experience an “**addictive**” effect due to which the impulse will be felt less and less: if necessary the intensity can be increased by one level to counter this effect.

Position of electrodes: form a square over the painful area as shown in illustration 1.

A02 • TENS endorphinic (Medical program)

This type of stimulation produces two types of effects according to how the electrodes are positioned: positioning the electrodes in the dorsal region, see photo 08 in the positions manual, promotes the endogenous production of morphine-like substances capable of raising the pain perception threshold; positioning the electrodes to form a square over the painful area

as shown in illustration 1 produces a vascularizing effect. Vascularization increases arterial flow and consequently aids the removal of algogenic substances and helps to restore normal physiological conditions.

Treatment duration 20 minutes in a single phase, daily frequency.

Do not position the electrodes close to inflamed areas.

Intensity adjusted for good solicitation of the part stimulated, the sensation must be similar to that of a massage.

A03 • TENS sequential (Medical program)

The action is very similar to that of the A01 program. In this program, the impulse frequency and amplitude are varied during stimulation in order to counter the effect of the user becoming accustomed to the stimulation (continuous increases of intensity will not be required). Select the intensity at the start of the program and maintain the same level until the end of the program. The user's perception of the current may almost become nonexistent during the program; if this occurs, do not increase the intensity and wait until the program has finished. The TENS program is working correctly.

Program duration: 20 minutes per phase. Intensity: slightly above the threshold of perception.

Position of the electrodes: form a square over the painful area as shown in illustration 1 or refer to the positions illustrated from photo 25 to photo 33.

A04 • TENS maximum value (Medical program)

Very short duration, 3 minutes. Blocks pain impulses peripherally creating a proper anaesthetising effect in the area treated. This type of stimulation is suitable for injuries or bruises when rapid action is required. The **intensity** selected is the maximum tolerable value (well in excess of conventional tens, and therefore with considerable contraction of the muscles surrounding the area treated). That is the reason why such stimulation is undoubtedly the least tolerated but is extremely effective. This type of stimulation is not recommended for particularly sensitive people and in any case the electrodes should not be positioned in sensitive areas such as the face and genitals or close to wounds.

Position of electrodes: form a square above the painful area as shown in illustration 1.

A05 • TENS Burst (Medical program)

This stimulation produces a training TENS effect by using the conventional TENS frequency. Commonly used for pain therapy. The action is similar to that of the endorphinic TENS program, and the treatment time is 15 minutes for each phase.

Position of electrodes: form a square above the painful area as shown in illustration 1.

A06 • Cervical pain (Medical program) (only on T-ONE EVO II)

Specific program for the treatment of pain in the neck area.

The **intensity** should be adjusted to a level between the thresholds of perception and pain: the maximum intensity level is the moment in which the muscles surrounding the area treated begin to contract; over this limit stimulation does not become more effective, just more irritating, so it is best to stop before that point. The first benefits can be seen after 10 to 12 treatments carried out on a daily basis; proceed with the treatment until the symptoms pass. Position of electrodes: photo 25.

A06 Cervical pain		Total duration: 33 minutes
Electrodes Positioning: cervical zone (photo 25).		
Phase 1	Phase 2	
Tens microimpulses 30 min.	Decontracting 3 min.	

A07 • Lumbal pain (Medical program) (only on T-ONE EVO II)

Specific program for the treatment of pain in the lumbar area. The intensity should be adjusted to a level between the thresholds of perception and pain: the maximum **intensity** level is the moment in which the muscles surrounding the area treated begin to contract; over this limit stimulation does not become more effective, just more irritating, so it is best to stop before that point. The first benefits can be seen after 15 to 20 treatments carried out on a daily basis; proceed with the treatment until the symptoms pass. Program duration: 40 minutes.
Position of electrodes: **see photo 27 in the manual of positions.**

A07 Lumbal pain		Total duration: 35 minutes
Electrodes Positioning: lumbar area (photo 27).		
Phase 1	Phase 2	
Tens microimpulses 30 min.	Tens Burst 5 min.	

A08 • Sciatalgy (Medical program) (only on T-ONE EVO II)

Specific program for the treatment of pain in the sciatic area. The intensity should be adjusted to a level between the thresholds of perception and pain: the maximum **intensity** level is the moment in which the muscles surrounding the area treated begin to contract; over this limit stimulation does not become more effective, just more irritating, so it is best to stop before that point. The first benefits can be seen after 15 to 20 treatments carried out on a daily basis; proceed with the treatment until the symptoms pass. Program duration: 30 minutes.
Position of electrodes: **see photo 28 in the manual of positions.**

A09 • Distortion/Contusion (Sprains / Bruises) (Medical program) (only on T-ONE EVO II)

The program develops its effectiveness after this type of injury by inhibiting pain locally, producing three selectively acting, differentiated impulses. The **intensity** should be adjusted to a level between the thresholds of perception and pain:
Number of treatments: until pain is lessened, on a daily basis (even 2/3 times a day).

A8 Sprains/Bruises		Total duration: 30 minutes
Electrodes Positioning: ankle (photo 32).		
Phase 1	Phase 2	
Tens rapid 15 min.	Haematomas impulse 15 min.	

A10 • Inflammatory (Medical program) (only on T-ONE EVO II)

Program recommended for inflammatory conditions. To be applied until the inflammatory state is lessened (10-15 applications, once a day; the daily treatments can be doubled if required). Identify the area to be treated and position the electrodes as shown in illustration 1. Adjust the **intensity** until a tingling feeling is produced in the area treated; avoid contracting the surrounding muscles.
Program duration: 30 minutes.

A11 • Vascularization (Medical program) (only on T-ONE EVO II)

Has a vascularizing effect on the area treated. Vascularization increases arterial flow and consequently aids the removal of algogenic substances and helps to restore normal physiological conditions. Do not position the electrodes close to inflamed areas. Daily application is recommended, the number of applications is not defined; the program can be used to reduce pain.
Stimulation **intensity** should be between the thresholds of perception and slight discomfort.

Program duration: 20 minutes. Position of electrodes: **see photo 25 and 33 in the manual of positions.**

A12 • Decontracting (only on T-ONE EVO II)

Can be used for all sports, after competitions or the most demanding training sessions, in particular after long and intense exertion. To be used immediately after exertion. Helps drainage and capillarization, improving muscle oxygenation and helping to discharge synthetic substances produced during exertion. Program duration: 35 minutes. Stimulation **intensity**: medium-low, increased during the last 10 minutes.

Position of electrodes: **see photos 01 to 20 in the manual of positions.**

A12 Decontracting		Total duration: 35 minutes	
Electrodes Positioning: form a square above the painful area as shown in illustration 1.			
Fase 1		Fase 2	
Tens a impulso stretto 15 min.		Tens a impulso largo 20 min.	

A13 • Hematoma (Medical program) (only on T-ONE EVO II)

Consult a doctor before using this program to treat haematomas. Total program duration: 20 minutes. Few applications carried out within a few hours of the bruise. A combination of various types of square-wave impulses has a graduated draining effect on the area to be treated (impulses at different frequencies drain the area at different depths). The **intensity** should be adjusted to a level between the thresholds of perception and pain, without causing muscle contraction.

Position of electrodes: form a square above the area to be treated as shown in illustration 1.

A14 • Regenerating skin (only on T-ONE EVO II)

This program has four phases with three different impulses to vascularize and tone the superficial muscles. It also produces a local pain-relieving action. The intensity should be adjusted to a level between the thresholds of perception and pain.

A14 Regenerating skin				Total duration: 40 minutes	
Electrodes Positioning: form a square above the area to treat as shown in illustration 1.					
Phase 1		Phase 2		Phase 3	
10 min. Vascularization		10 min. Tens micropulses		10 min. Tens sequential	
				Phase 4	
				10 min. Vascularization	

A15 • Atrophy prevention (Medical program) (only on T-ONE EVO II)

Program created to maintain muscle trophism.

This treatment concentrates on muscle toning, paying particular attention to slow-twitch fibres. Particularly indicated for patients recovering from an accident or an operation. Prevents muscle atrophy caused by physical inactivity. The muscle area concerned can be stimulated with daily applications of medium **intensity**; if you increase the intensity, leave a day of rest between applications to allow the muscles to recover. The intensity of this program must be adjusted to produce, during second phase, good muscle contraction in the area treated.

Position of electrodes from photo 1 to photo 20.

Program duration: 30 minutes.

A15 Atrophy prevention		Total duration: 30 minutes	
			
Biceps (photo 02/15), Triceps (photo 03/16), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).			

	Abdominals (photo 01/20), Pectoral/Breast (photo 07/17), Trapezius (photo 08), Large dorsal (photo 09), Gluteus (photo 19)		
	Quadriceps/Thighs (photo 11/18), Biceps femoris (photo12), Calves (photo13), Anterior Tibial (photo14).		
Phase 1	Phase 2	Phase 3	Phase 4
5 min. warm up	10 min. 9 sec recovery 6 sec contraction	10 min. 9 sec recovery 6 sec contraction	5 min. recovery

A16 • Atrophy (trophism rehabilitation) (Medical program) (only on T-ONE EVO

II)

This program acts selectively on slow-twitch fibres. Ideal for recovering muscle trophism after a long period of inactivity or an accident.

Program to be carried out when loss of muscle tone has already occurred. Apply with caution (at low **intensity**, enough to produce light muscle contractions) in the first 2/3 weeks. Increase intensity progressively over the next 3/4 weeks. Application on alternate days. Position of electrodes from photo 1 to photo 20.

Program duration: 30 minutes.

A16 Atrophy		Total duration: 30 minutes
	Biceps (photo 02/15), Triceps (photo 03/16), hand extensors (photo 04), hand flexors (photo 05), Deltoid (photo 06).	
	Abdominals (photo 01/20), Pectoral/Breast (photo 07/17), Trapezius (photo 08), Large dorsal (photo 09), Gluteus (photo 19)	
	Quadriceps/Thighs (photo 11/18), Biceps femoris (photo12), Calves (photo13), Anterior Tibial (photo14).	
Phase 1	Phase 2	Phase 3
3 min. warm up	25 min. 6 sec. Recovery 6 sec. contraction	2 min. recovery

14. SET UP.

Selecting the language and contrast level on the display.

Press the ON/OFF button  and hold it down until the language list appears on the display. Release the buttons: now use the mouse to move the cursors and choose the desired language or the level of contrast on the display.

Select the chosen language by using the   buttons.

Adjust the level of contrast on the display by using the   buttons.

Press  to confirm your selections.

15. Self-adhesive electrodes.

The position of the electrodes is shown in the manual of electrode positions.

The electrodes supplied are top quality, pre-gelled and ready for use, they are particularly suitable for electrostimulation treatments.

They are very flexible making it easy to treat all areas. Remove the protective plastic from the electrode before use, position it on the skin as shown in the manual and replace the plastic after use.

Electrode duration is determined by the skin's pH value; we recommend no more than 20 applications using the same electrode.

Repeated use of the same electrodes can compromise stimulation safety: do not use electrodes when they no longer stick to the skin; in fact it can cause skin redness that can last for many hours after stimulation. Consult a doctor if this occurs.


16. Declaration of Conformity

IACER S.r.l., headquartered in Italy, via S. Pertini 24/A 30030 Martellago (VE), declares on its own responsibility devices T-ONE EVO I and T-ONE EVO II is manufactured in conformity with 93/42/EEC Directive (MDD) as amended by 2007/47/CE date 25 September 2007.

T-ONE EVO I and T-ONE EVO II are Class IIa devices, with reference to 93/42/EEC Directive (MDD), annex IX rule 9 (and following modifications Dir. 2007/47/CE)

Martellago, 29/03/2010

Legal representative
Mario Caprara

17. Technical specifications.

Risk class: IIa in accordance with the 93/42/EEC directive (and following modifications).

Insulation class: Internally powered device with type BF applied parts, in accordance with the EN 60601-1 standard.

Classification: IP30, based on the penetration of liquids and powder; equipment not protected.

Housing: ABS.

Power Supply: Ni/Mh rechargeable battery 4.8V 1600mAh.

Battery charger: model XKD-C0300NHS6.0-5W, input 100-240 VAC 50-60Hz, output 6 VDC 300mA or other compatible model supply by manufacturer

Ave. consumption: 1 W.

Limitations: equipment not suitable for use in environments with a high concentration of oxygen and/or flammable liquids and/or flammable gas; do not use with equipment for electrosurgery or short-wave or microwave therapy.

Operation: equipment suitable for continuous operation.

Output channel: 4 channels, independent and galvanically isolated.

Output current: maximum power, peak to peak 240mA, 1Kohm load.

Type of impulse: compensated biphasic square wave.

Frequency: operates with a range of impulses from 2 Hz to 660 Hz (Tens Burst).

Impulse amplitude: from 60µs to 350µs.

Display: reflective LCD display.

Control keypad: 10-key keypad in ABS.

Dimensions: 125X76X34 mm.

Weight: 280 gr. including batteries.

Maintenance and storage: refer to the instructions in Chapter 8 "Maintenance and storage of T-ONE".

Warning: the equipment delivers current in excess of 10mA.

CE 0476 Produced in compliance with the 93/42/EEC directive concerning medical devices (and following modifications Dir. 2007/47/CE).

18. Assistance

Every intervention on device must be performed by manufacturer. For any assistance intervention contact the National Distributor or the manufacturer at the following address:

I.A.C.E.R. S.r.l.

Via S. Pertini, 24/a • 30030 Martellago (VE) - Italy

Tel. +39 0415401356 • Fax +39 0415402684

www.iacer.it • iacer@iacer.it

You can get any technical documentation on spare parts but only prior business authorization.

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I.A.C.E.R S.r.l. SU